

β -actin(5B7)Mouse Monoclonal Antibody

产品货号	GP160-100
产品名称	β -actin(5B7)Mouse Monoclonal Antibody
类别	抗原抗体
推荐应用	WB, IHC-p, Elisa
反应种属	Human, Rat, Mouse, Mk, Dg, Ch, Hamster, Rb, Insect
浓度	1mg/ml
存储缓冲液	Liquid in PBS,pH 7.4,containing 50%glycerol,2%BSA and New type preservative N
免疫原	Synthetic Peptide of β -actin
特异性	The antibody detects endogenous β -actin protein.
稀释度	IF:1:50-200 WB:1:5000 IHC:1:200
参考分子量	42KD
宿主	Mouse
纯化	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen