

### GAPDH Rabbit Polyclonal Antibody

产品货号	GP150-100
产品名称	GAPDH Rabbit Polyclonal Antibody
类别	抗原抗体
推荐应用	WB, IHC-p, Elisa
反应种属	Human, Mouse, Rat, Rb, Ch, Mk, sheep, Xenopus
浓度	1mg/ml
存储缓冲液	Liquid in PBS, pH 7.4, containing 50% glycerol, 2% BSA and New type preservative N.
免疫原	Recombinant Protein of GAPDH
特异性	The antibody detects endogenous GAPDH protein
稀释度	WB:1:5000 IHC:1:200
参考分子量	37KD
宿主	Rabbit
背景介绍	Glyceraldehyde 3 phosphate dehydrogenase (GAPDH) is well known as one of the key enzymes involved in glycolysis. GAPDH is constitutively expressed in almost all tissues at high levels, therefore antibodies against GAPDH are useful as loading controls
纯化	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using specific immunogen